

NeuroTrace Study Guide

Domain: Domain I – Basic Concepts & Principles

Section: EEG Electrode Placement: The 10–20 System

Style: Spatial, applied, exam-oriented

1. Core Principles (Must Know)

The 10–20 System Is Based on Relative Skull Measurements

- **The 10–20 system is based on relative skull measurements**
- Uses percentages of skull dimensions
- Not absolute distances
- Adapts to different head sizes
- Ensures proportional placement

Ensures Reproducible Electrode Placement

- **Ensures reproducible electrode placement**
- Standardized system
- Consistent across patients
- Reproducible across studies
- Enables comparison

Supports Consistent Localization Across Patients

- **Supports consistent localization across patients**
- Same electrode = same brain region
- Enables comparison between patients
- Supports research and clinical practice
- Essential for interpretation

Key Principle

- **Standardized placement enables meaningful EEG interpretation**
- Without standardization, interpretation is unreliable
- Standardized placement is essential
- Enables accurate localization
- Supports clinical correlation

Practical Application

- Use 10–20 system for all recordings
 - Follow measurement procedures
 - Verify electrode placement
 - Document placement accuracy
-

2. Naming Conventions

Letter = Brain Region

F = Frontal

- **F = Frontal:** Frontal lobe electrodes
- Fp1, Fp2: Frontal poles

- F3, F4: Frontal regions
- F7, F8: Anterior temporal/frontal
- Fz: Frontal midline

T = Temporal

- **T = Temporal:** Temporal lobe electrodes
- T3, T4: Mid-temporal
- T5, T6: Posterior temporal
- T7, T8: Alternative naming (same as T3, T4)
- Critical for epilepsy localization

C = Central

- **C = Central:** Central (Rolandic) region
- C3, C4: Central regions
- Cz: Central midline
- Over motor/sensory cortex
- Mu rhythm location

P = Parietal

- **P = Parietal:** Parietal lobe electrodes
- P3, P4: Parietal regions
- Pz: Parietal midline
- Less commonly used
- May be used in extended montages

O = Occipital

- **O = Occipital:** Occipital lobe electrodes
- O1, O2: Occipital regions
- Oz: Occipital midline
- Alpha rhythm (PDR) location
- Visual cortex region

Numbers

Odd = Left Hemisphere

- **Odd = Left hemisphere:** 1, 3, 5, 7
- Fp1, F3, C3, P3, O1: Left side
- F7, T3, T5: Left temporal chain
- Critical for localization

Even = Right Hemisphere

- **Even = Right hemisphere:** 2, 4, 6, 8
- Fp2, F4, C4, P4, O2: Right side
- F8, T4, T6: Right temporal chain
- Critical for localization

"z" = Midline

- **"z" = Midline:** Central placement
- Fz: Frontal midline
- Cz: Central midline
- Pz: Parietal midline
- Oz: Occipital midline

Best Practice

- Memorize letter meanings (F, T, C, P, O)
 - Remember odd = left, even = right
 - Know that "z" = midline
 - Use naming to identify locations
-

3. Electrode Locations & Regions

Frontal Poles: Fp1, Fp2

- **Frontal poles: Fp1, Fp2**
- Most anterior electrodes
- Over frontal poles
- May show eye movement artifacts
- Important for frontal lobe assessment

Central Midline: Cz

- **Central midline: Cz**
- Top of head (vertex)
- Over central sulcus
- Mu rhythm location
- Reference point for measurements

Posterior Dominant Regions: O1, O2

- **Posterior dominant regions: O1, O2**
- Over occipital cortex
- Alpha rhythm (PDR) location
- Visual cortex region
- Critical for background assessment

Temporal Chains: F7–T3–T5 / F8–T4–T6

- **Temporal chains: F7–T3–T5 / F8–T4–T6**
- F7/F8: Anterior temporal/frontal
- T3/T4: Mid-temporal
- T5/T6: Posterior temporal
- Critical for epilepsy localization

ABRET Emphasis

- **Temporal electrodes are critical for epilepsy localization**
- Most epileptiform activity is temporal
- Temporal electrodes are essential
- Must know temporal chain locations
- Critical for focal epilepsy

Best Practice

- Know electrode locations
 - Understand regional groupings
 - Recognize temporal chain importance
 - Use locations for localization
-

4. Measurement Logic

Distances Measured as Percentages

10% and 20% Intervals

- **10% and 20% intervals:** Measurement basis
- 10%: Between adjacent electrodes
- 20%: Between major electrodes
- Based on skull dimensions
- Ensures proportional placement

Reference Points

Nasion

- **Nasion:** Bridge of nose
- Anterior reference point
- Used for frontal measurements
- Landmark for placement

Inion

- **Inion:** Back of head (occipital protuberance)
- Posterior reference point
- Used for occipital measurements
- Landmark for placement

Preauricular Points

- **Preauricular points:** In front of ears
- Lateral reference points
- Used for temporal measurements
- Landmark for placement

Best Practice

- Understand measurement logic
- Know reference points
- Use percentages, not absolute distances
- Verify measurements

5. Clinical Localization Examples

EEG Finding	Likely Region	Electrodes
Focal spikes at T3	Left temporal	T3, F7, T5
Slowing at F4	Right frontal	F4, Fp2, F8
PDR maximal at O1/O2	Occipital	O1, O2, Oz
Spikes at C3	Left central	C3, F3, P3
Generalized spike-and-wave	Both hemispheres	All electrodes

Key Distinctions

Focal Findings

- **Focal findings:** Limited to specific region
- One or few electrodes involved
- Localized to specific area
- Suggests focal pathology

Generalized Findings

- **Generalized findings:** Both hemispheres
- Multiple electrodes involved
- Widespread distribution
- Suggests generalized process

Best Practice

- Use electrode names for localization
 - Know which electrodes = which regions
 - Understand focal vs generalized
 - Apply localization logic
-

6. Common ABRET Exam Traps

Trap 1: Confusing Left vs Right Numbering

- **Reality:** Odd = left, even = right
- **Trap:** May confuse left and right
- **Solution:** Memorize odd = left, even = right
- **ABRET focus:** Electrode naming

Trap 2: Misidentifying Electrode Regions

- **Reality:** Letters indicate regions
- **Trap:** May misidentify regions
- **Solution:** Memorize letter meanings
- **ABRET focus:** Regional identification

Trap 3: Assuming Cz Is Frontal

- **Reality:** Cz is central, not frontal
- **Trap:** May think Cz is frontal
- **Solution:** Know that C = central
- **ABRET focus:** Electrode locations

Trap 4: Forgetting Midline Designation

- **Reality:** "z" = midline
- **Trap:** May forget midline designation
- **Solution:** Remember "z" = midline
- **ABRET focus:** Midline electrodes

Trap 5: Not Understanding Temporal Chain

- **Reality:** Temporal chain is F7–T3–T5 / F8–T4–T6
 - **Trap:** May not understand temporal chain
 - **Solution:** Learn temporal chain sequence
 - **ABRET focus:** Temporal localization
-

7. Clinical Correlation

Accurate Placement

Improves Localization

- **Improves localization:** Accurate placement = accurate localization
- Correct electrode position = correct brain region
- Enables precise localization
- Supports clinical correlation

Reduces Misinterpretation

- **Reduces misinterpretation:** Accurate placement prevents errors
- Misplacement leads to false localization
- Accurate placement prevents confusion
- Essential for correct interpretation

Poor Placement Leads to False Localization

- **Poor placement leads to false localization**
- Incorrect electrode position = incorrect localization
- May lead to wrong diagnosis
- May lead to wrong treatment
- Must ensure accurate placement

Best Practice

- Ensure accurate electrode placement
 - Verify measurements
 - Document placement accuracy
 - Correct misplacement immediately
-

8. Case-Based Example

Scenario

Clinical Setting: Routine EEG for seizure evaluation

EEG Finding: Spike maximal at F7

Question: What region is this?

Interpretation

- **Left anterior temporal/frontal region**
- F7 = Left anterior temporal/frontal
- Spike maximal at F7 = focal finding
- Suggests left anterior temporal/frontal pathology
- May indicate focal epilepsy

Teaching Point

- **Electrode name encodes both side and region**
- F7 tells us: F = frontal/temporal, 7 = left, anterior
- Electrode name provides localization
- Essential for interpretation
- Must know naming conventions

ABRET Application

- Given electrode name → identify side and region
 - Given F7 → left anterior temporal/frontal
 - Given T3 → left mid-temporal
 - Must know naming conventions
-

9. Exam Readiness Checklist

Use this checklist to verify your understanding:

- ☐ Can identify electrode locations (frontal, temporal, central, parietal, occipital)
 - ☐ Can apply numbering rules (odd = left, even = right)
 - ☐ Can use placement for localization (electrode name → brain region)
 - ☐ Can avoid left/right confusion (odd vs even)
 - ☐ Understand that 10–20 system is based on relative measurements
 - ☐ Know that standardized placement enables interpretation
 - ☐ Recognize that temporal electrodes are critical for epilepsy
 - ☐ Know that "z" = midline
 - ☐ Can identify electrode regions from names
 - ☐ Understand measurement logic (10%, 20%, reference points)
 - ☐ Can localize findings based on electrode positions
 - ☐ Can identify common ABRET exam traps
-

10. Internal Cross-Links

Foundations

- **Neuroanatomy for EEG Localization:** Understanding brain regions
- **Normal EEG Rhythms:** Where rhythms are recorded

Workflow

- **Montages & Referencing:** How electrodes are connected
- **Electrodes & Impedance:** Electrode application and maintenance

Patterns

- **Epileptiform Discharges:** Localization of epileptiform activity
- **Focal vs Generalized:** Using electrodes for classification

Cases

- **Localization-based EEG cases:** Cases requiring electrode knowledge
- **Montage interpretation cases:** Cases using electrode positions

Quizzes

- **10–20 system MCQs:** Questions on electrode placement
 - **Localization questions:** Questions on using electrodes for localization
 - **Naming convention questions:** Questions on electrode names
-

Study Tips

1. **Memorize naming conventions:** F, T, C, P, O = regions; odd = left, even = right; "z" = midline
 2. **Learn electrode locations:** Frontal poles, temporal chains, central, occipital
 3. **Know measurement logic:** 10%, 20% intervals, reference points
 4. **Understand localization:** Electrode name → brain region
 5. **Remember the principle:** Standardized placement enables interpretation
 6. **Know the traps:** Left/right confusion, regional misidentification, midline forgetting
 7. **ABRET focus:** Expect questions on naming, locations, and localization
-

End of Study Guide

For additional practice, complete quiz questions tagged: electrode-placement, 10-20, localization, hemispheres