

# NeuroTrace Academy Study Guide

**Domain:** Domain II – Performing the EEG Study

**Section:** Complete Electrode Placement Guide: Standard & Exotic Electrodes

**Style:** Comprehensive, measurement-based, clinical application, exam-oriented

---

## 1. Core Principles (Must Know)

### Standard 10-20 System Is the Foundation

- **10-20 system:** International standard for routine EEG
- **21 standard electrodes:** Fp1, Fp2, F3, F4, F7, F8, Fz, C3, C4, Cz, P3, P4, Pz, T3, T4, T5, T6, O1, O2, Oz, plus ground/reference
- **Based on percentages:** 10% and 20% of head measurements
- **Reproducible:** Consistent placement across patients and studies

### Extended Systems for Specialized Applications

- **10-10 system:** Extended system with more electrodes (additional 10% positions)
- **Exotic electrodes:** Additional electrodes for specific clinical needs
- **Invasive electrodes:** For presurgical evaluation (sphenoidal, nasopharyngeal, depth, subdural)

### Key Principle

#### Electrode selection depends on clinical indication

- Routine EEG: Standard 10-20 system
- Temporal lobe epilepsy: Add T1, T2, or sphenoidal electrodes
- Presurgical evaluation: May require invasive electrodes
- Research: May use extended 10-10 system

### Practical Application

- Know standard 10-20 system thoroughly
  - Understand when to add exotic electrodes
  - Know placement techniques for all electrode types
  - Understand clinical indications for each
- 

## 2. Standard 10-20 System (Review)

# An Introduction to Brain Structures




Figure 1: Lateral view of the brain showing the relationship between brain anatomy and EEG electrode placement. The precentral gyrus (motor cortex, red) corresponds to C3/C4 electrodes. The postcentral gyrus (sensory cortex, blue) also corresponds to C3/C4. Frontal electrodes (Fp1, F3, F7) overlie the frontal lobe, temporal electrodes (T3, T5) overlie the temporal lobe, and occipital electrodes (O1, O2) overlie the occipital lobe.

## Standard Electrode Set (21 Electrodes)

### Frontal:

- Fp1, Fp2: Frontal poles (10% above nasion)
- F3, F4: Frontal regions (20% from midline)
- F7, F8: Anterior temporal/frontal (10% above Fp, at temporal line)
- Fz: Frontal midline

### Central:

- C3, C4: Central regions (20% from midline)

- Cz: Central midline (vertex)

#### **Parietal:**

- P3, P4: Parietal regions (20% from midline)
- Pz: Parietal midline

#### **Temporal:**

- T3, T4: Mid-temporal (at temporal line)
- T5, T6: Posterior temporal

#### **Occipital:**

- O1, O2: Occipital regions (10% above inion, 20% from midline)
- Oz: Occipital midline

#### **Reference:**

- A1: Left ear (mastoid or earlobe)
- A2: Right ear (mastoid or earlobe)

#### **Ground:**

- G: Ground electrode (typically Fpz or separate site)

### **Naming Conventions**

- **Letters:** F=Frontal, T=Temporal, C=Central, P=Parietal, O=Occipital
- **Numbers:** Odd=Left, Even=Right
- **"z":** Midline electrodes (Fz, Cz, Pz, Oz)

### **Measurement Reference Points**

1. **Nasion:** Indentation at top of nose (between eyebrows)
2. **Inion:** Most prominent point of occipital bone
3. **Preauricular points:** Indentations anterior to tragus of each ear
4. **Cz (Vertex):** Midpoint of nasion-inion and preauricular lines

## **3. Extended 10-10 System**

### **Additional Electrode Positions**

**10-10 system adds electrodes at 10% intervals between standard 20% positions:**

#### **Additional Frontal:**

- AF3, AF4, AF7, AF8, AFz: Anterior frontal positions
- FC1, FC2, FC3, FC4, FC5, FC6, FCz: Frontal-central positions

#### **Additional Central:**

- C1, C2, C5, C6: Additional central positions

#### **Additional Parietal:**

- CP1, CP2, CP3, CP4, CP5, CP6, CPz: Central-parietal positions
- P1, P2, P5, P6, P9, P10: Additional parietal positions

#### **Additional Temporal:**

- FT7, FT8, FT9, FT10: Frontal-temporal positions
- TP7, TP8, TP9, TP10: Temporal-parietal positions

**Additional Occipital:**

- PO3, PO4, PO7, PO8, POz: Parietal-occipital positions

**Clinical Use**

- **Research:** High-density EEG studies
  - **Source localization:** More electrodes improve spatial resolution
  - **Presurgical evaluation:** Better localization of epileptogenic zones
  - **Not routine:** Typically used in specialized settings
- 

## 4. Exotic Electrodes: T1 and T2

### T1 and T2 Electrodes

**Location:**

- **T1:** Left inferior temporal (1 cm above left ear canal, 1/3 distance from T3 to external auditory meatus)
- **T2:** Right inferior temporal (1 cm above right ear canal, 1/3 distance from T4 to external auditory meatus)

**Placement Technique:**

1. Locate T3 (T7) or T4 (T8) position
2. Locate external auditory meatus (ear canal)
3. Place electrode 1/3 of distance from T3/T4 toward ear canal
4. Approximately 1 cm above ear canal

**Clinical Indication:**

- **Temporal lobe epilepsy:** Better detection of mesial temporal activity
- **Presurgical evaluation:** Improved localization of temporal foci
- **When standard T3/T4 miss activity:** Add T1/T2 for better coverage

**Advantages:**

- Closer to temporal lobe than T3/T4
- Better detection of mesial temporal activity
- Non-invasive (surface electrode)
- Easy to apply

**Disadvantages:**

- May pick up more muscle artifact (near temporalis muscle)
- May pick up more EKG artifact
- Not as sensitive as sphenoidal electrodes

**ABRET Exam Focus:**

- T1/T2 are inferior temporal electrodes
  - Used for temporal lobe epilepsy
  - Placement: 1/3 distance from T3/T4 to ear canal
- 

## 5. Exotic Electrodes: F9 and F10

## F9 and F10 Electrodes

### Location:

- **F9:** Left inferior frontal (below F7, near left temple)
- **F10:** Right inferior frontal (below F8, near right temple)

### Placement Technique:

1. Locate F7 (left) or F8 (right) position
2. Place electrode inferior (below) F7/F8
3. Approximately at level of zygomatic arch
4. Near temple region

### Clinical Indication:

- **Frontal lobe epilepsy:** Better detection of inferior frontal activity
- **Orbitofrontal seizures:** Improved coverage of orbitofrontal cortex
- **When standard F7/F8 miss activity:** Add F9/F10 for better coverage

### Advantages:

- Closer to inferior frontal cortex
- Better detection of orbitofrontal activity
- Non-invasive (surface electrode)

### Disadvantages:

- May pick up more muscle artifact (temporalis, masseter)
- May pick up more eye movement artifact
- Less commonly used than T1/T2

### ABRET Exam Focus:

- F9/F10 are inferior frontal electrodes
  - Used for frontal lobe epilepsy, especially orbitofrontal
  - Placement: Below F7/F8, near temple
- 

## 6. Exotic Electrodes: P9 and P10

### P9 and P10 Electrodes

#### Location:

- **P9:** Left inferior parietal (posterior to T5, near left mastoid)
- **P10:** Right inferior parietal (posterior to T6, near right mastoid)

#### Placement Technique:

1. Locate T5 (left) or T6 (right) position
2. Place electrode posterior and inferior to T5/T6
3. Near mastoid region
4. Posterior temporal/inferior parietal location

#### Clinical Indication:

- **Posterior temporal/parietal epilepsy:** Better detection of inferior parietal activity
- **Presurgical evaluation:** Improved coverage of posterior regions

- **Less commonly used:** Not as common as T1/T2 or F9/F10

**Advantages:**

- Better coverage of inferior parietal cortex
- Non-invasive (surface electrode)

**Disadvantages:**

- Less commonly needed
- May pick up more muscle artifact
- Less sensitive than specialized electrodes

**ABRET Exam Focus:**

- P9/P10 are inferior parietal electrodes
  - Used for posterior temporal/parietal epilepsy
  - Placement: Posterior and inferior to T5/T6
- 

## 7. Sphenoidal Electrodes (Sp1, Sp2)

### Sphenoidal Electrodes

**Location:**

- **Sp1:** Left sphenoidal (inserted through cheek, near foramen ovale)
- **Sp2:** Right sphenoidal (inserted through cheek, near foramen ovale)

**Placement Technique:**

1. **Insertion site:** Through cheek, anterior to mandibular condyle
2. **Depth:** Advanced to near foramen ovale (medial to mandible)
3. **Method:** Needle electrode inserted percutaneously
4. **Position:** Near mesial temporal structures

**Clinical Indication:**

- **Temporal lobe epilepsy:** Gold standard for mesial temporal detection
- **Presurgical evaluation:** Essential for temporal lobe epilepsy surgery
- **When surface electrodes miss activity:** Sphenoidal electrodes are more sensitive

**Advantages:**

- **Most sensitive for mesial temporal activity**
- Closer to hippocampus and amygdala
- Better detection of temporal lobe seizures
- Gold standard for temporal lobe epilepsy

**Disadvantages:**

- **Invasive:** Requires needle insertion
- **Discomfort:** Patient may experience pain or discomfort
- **Risk:** Small risk of infection, bleeding
- **Not routine:** Only for specialized evaluations

**Placement Details:**

- Inserted through cheek (percutaneous)
- Advanced to near foramen ovale

- Positioned medial to mandible
- Near mesial temporal structures

**ABRET Exam Focus:**

- Sphenoidal electrodes are invasive
  - Used for temporal lobe epilepsy (presurgical evaluation)
  - Most sensitive for mesial temporal activity
  - Inserted through cheek, near foramen ovale
- 

## 8. Nasopharyngeal Electrodes (NP1, NP2)

### Nasopharyngeal Electrodes

**Location:**

- **NP1:** Left nasopharyngeal (inserted through nostril, advanced to nasopharynx)
- **NP2:** Right nasopharyngeal (inserted through nostril, advanced to nasopharynx)

**Placement Technique:**

1. **Insertion:** Through nostril, advanced to nasopharynx
2. **Depth:** Advanced until tip is in nasopharynx (behind soft palate)
3. **Method:** Flexible wire electrode with tip
4. **Position:** Near mesial temporal structures (via nasopharynx)

**Clinical Indication:**

- **Temporal lobe epilepsy:** Alternative to sphenoidal electrodes
- **Presurgical evaluation:** When sphenoidal not available or contraindicated
- **Less invasive than sphenoidal:** But less sensitive

**Advantages:**

- Less invasive than sphenoidal
- Better than surface electrodes for mesial temporal activity
- Can be used when sphenoidal contraindicated

**Disadvantages:**

- **Less sensitive than sphenoidal:** Not as close to mesial temporal structures
- **Discomfort:** Patient may experience nasal discomfort
- **Artifact:** May pick up swallowing, respiratory artifact
- **Less commonly used:** Sphenoidal preferred when possible

**Placement Details:**

- Inserted through nostril
- Advanced to nasopharynx (behind soft palate)
- Flexible wire electrode
- Near mesial temporal structures

**ABRET Exam Focus:**

- Nasopharyngeal electrodes are less invasive than sphenoidal
- Used for temporal lobe epilepsy
- Inserted through nostril to nasopharynx

- Less sensitive than sphenoidal but better than surface
- 

## 9. Additional Specialized Electrodes

### EOG (Electrooculogram) Electrodes

#### Location:

- **EOG1:** Above left eye (for vertical eye movements)
- **EOG2:** Below left eye (for vertical eye movements)
- **EOG3:** Lateral to left eye (for horizontal eye movements)
- **EOG4:** Lateral to right eye (for horizontal eye movements)

#### Purpose:

- Detect eye movements
- Distinguish eye movement artifact from brain activity
- Monitor eye movements during sleep studies

### EMG (Electromyogram) Electrodes

#### Location:

- **Chin EMG:** Under chin (for muscle tone during sleep)
- **Limb EMG:** On limbs (for movement detection)

#### Purpose:

- Monitor muscle activity
- Detect movements during sleep
- Distinguish muscle artifact from brain activity

### EKG (Electrocardiogram) Electrodes

#### Location:

- **EKG:** On chest or limb (for heart rate monitoring)

#### Purpose:

- Monitor heart rate
  - Detect EKG artifact in EEG
  - Essential for long-term monitoring
- 

## 10. Clinical Indications Summary

### When to Use Standard 10-20 System

- **Routine EEG:** Standard evaluation
- **Initial seizure evaluation:** First-line assessment
- **Follow-up studies:** Comparison with previous studies
- **Most clinical indications:** Standard system is sufficient

### When to Add T1/T2

- **Temporal lobe epilepsy:** Suspected or confirmed
- **Mesial temporal activity suspected:** When T3/T4 miss activity



- **Presurgical evaluation:** Before temporal lobe surgery
- **Refractory temporal seizures:** When standard electrodes insufficient

### When to Add F9/F10

- **Frontal lobe epilepsy:** Suspected or confirmed
- **Orbitofrontal seizures:** Suspected inferior frontal activity
- **When F7/F8 miss activity:** Add for better coverage

### When to Use Sphenoidal Electrodes

- **Temporal lobe epilepsy (presurgical):** Essential for surgery planning
- **Refractory temporal seizures:** When surface electrodes insufficient
- **Mesial temporal sclerosis:** Suspected or confirmed
- **Presurgical evaluation:** Before temporal lobe resection

### When to Use Nasopharyngeal Electrodes

- **Temporal lobe epilepsy:** When sphenoidal contraindicated
- **Alternative to sphenoidal:** When sphenoidal not available
- **Less invasive option:** For patients who cannot tolerate sphenoidal

### When to Use Extended 10-10 System

- **Research studies:** High-density EEG
- **Source localization:** Precise localization needed
- **Presurgical evaluation:** Advanced cases
- **Not routine:** Specialized applications only

## 11. Placement Techniques

### Standard 10-20 System Placement

1. **Measure head:**
  - Nasion to inion (sagittal)
  - Left to right preauricular (coronal)
  - Head circumference
2. **Calculate positions:**
  - Use 10% and 20% increments
  - Mark electrode positions
3. **Apply electrodes:**
  - Prepare skin (abrade if needed)
  - Apply conductive paste/gel
  - Secure electrodes
  - Check impedance (<5 kΩ, ideally <2 kΩ)

### T1/T2 Placement

1. **Locate T3 (T7) or T4 (T8):** Standard mid-temporal position
2. **Locate external auditory meatus:** Ear canal opening
3. **Measure distance:** From T3/T4 to ear canal
4. **Place electrode:** 1/3 distance from T3/T4 toward ear canal

5. **Verify:** Approximately 1 cm above ear canal

### F9/F10 Placement

1. **Locate F7 (left) or F8 (right):** Standard anterior temporal/frontal position
2. **Place electrode:** Inferior (below) F7/F8
3. **Position:** Near zygomatic arch, temple region
4. **Verify:** Below F7/F8, near temple

### P9/P10 Placement

1. **Locate T5 (left) or T6 (right):** Standard posterior temporal position
2. **Place electrode:** Posterior and inferior to T5/T6
3. **Position:** Near mastoid region
4. **Verify:** Posterior temporal/inferior parietal location

### Sphenoidal Electrode Placement

1. **Identify insertion site:** Through cheek, anterior to mandibular condyle
2. **Insert needle:** Percutaneous insertion through cheek
3. **Advance electrode:** To near foramen ovale (medial to mandible)
4. **Verify position:** Near mesial temporal structures
5. **Secure:** Tape or secure externally
6. **Monitor:** Check for complications (bleeding, infection)

**Note:** Sphenoidal placement requires specialized training and is typically performed by physicians or specially trained technologists.

### Nasopharyngeal Electrode Placement

1. **Prepare electrode:** Flexible wire with tip
2. **Insert through nostril:** Gently advance through nasal passage
3. **Advance to nasopharynx:** Behind soft palate
4. **Verify position:** Tip in nasopharynx
5. **Secure:** Tape to nose/face externally
6. **Monitor:** Check for discomfort, displacement

---

## 12. Artifact Considerations

### T1/T2 Artifacts

- **Muscle artifact:** Temporalis muscle (near electrode)
- **EKG artifact:** Heart activity may be detected
- **Movement artifact:** Jaw movements

#### Mitigation:

- Proper placement (away from muscle)
- Patient instruction (relax jaw)
- Filtering if necessary

### F9/F10 Artifacts

- **Muscle artifact:** Temporalis, masseter muscles
- **Eye movement artifact:** Blinks, eye movements
- **Jaw movement artifact:** Talking, chewing

#### Mitigation:

- Proper placement
- Patient instruction
- Filtering if necessary

### Sphenoidal Artifacts

- **Muscle artifact:** Less than surface electrodes
- **Movement artifact:** Jaw movements may affect
- **Swallowing artifact:** May be detected

#### Mitigation:

- Proper placement
- Patient instruction
- Less artifact than surface electrodes

### Nasopharyngeal Artifacts

- **Swallowing artifact:** Common
- **Respiratory artifact:** Breathing may affect
- **Movement artifact:** Head movements

#### Mitigation:

- Proper placement
- Patient instruction
- Secure placement

## 13. ABRET Exam High-Yield Topics

### Must-Know Facts

1. **T1/T2:** Inferior temporal electrodes, 1/3 distance from T3/T4 to ear canal
2. **F9/F10:** Inferior frontal electrodes, below F7/F8
3. **P9/P10:** Inferior parietal electrodes, posterior to T5/T6
4. **Sphenoidal:** Invasive, through cheek, near foramen ovale, most sensitive for mesial temporal
5. **Nasopharyngeal:** Less invasive, through nostril, less sensitive than sphenoidal
6. **Clinical indications:** Know when to use each electrode type

### Common Exam Questions

#### Question Type 1: Electrode Location

- "T1 electrode is placed..."
- Answer: 1/3 distance from T3 to ear canal, inferior temporal

#### Question Type 2: Clinical Indication

- "Sphenoidal electrodes are used for..."
- Answer: Temporal lobe epilepsy, presurgical evaluation, mesial temporal activity

#### Question Type 3: Sensitivity

- "Which electrode is most sensitive for mesial temporal activity?"
- Answer: Sphenoidal electrodes

#### Question Type 4: Invasiveness

- "Which is less invasive: sphenoidal or nasopharyngeal?"

- Answer: Nasopharyngeal (but less sensitive)

#### Question Type 5: Placement Technique

- "F9 electrode is placed..."
- Answer: Below F7, near temple, inferior frontal

## 14. Comparison Table

Electrode	Location	Invasiveness	Sensitivity (Temporal)	Clinical Use
<b>T3/T4</b>	Mid-temporal	Non-invasive	Low-Moderate	Standard, routine
<b>T1/T2</b>	Inferior temporal	Non-invasive	Moderate	Temporal epilepsy
<b>F7/F8</b>	Anterior temporal/frontal	Non-invasive	Low-Moderate	Standard, routine
<b>F9/F10</b>	Inferior frontal	Non-invasive	Moderate	Frontal epilepsy
<b>Sphenoidal</b>	Near foramen ovale	Invasive	<b>Highest</b>	Presurgical, temporal
<b>Nasopharyngeal</b>	Nasopharynx	Less invasive	Moderate-High	Alternative to sphenoidal

## 15. Exam Readiness Checklist

Use this checklist to verify your understanding:

- ☐ Know standard 10-20 system (all 21 electrodes)
- ☐ Understand extended 10-10 system (additional positions)
- ☐ Know T1/T2 placement (1/3 distance from T3/T4 to ear canal)
- ☐ Know F9/F10 placement (below F7/F8, near temple)
- ☐ Know P9/P10 placement (posterior to T5/T6)
- ☐ Understand sphenoidal electrodes (invasive, through cheek, near foramen ovale)
- ☐ Understand nasopharyngeal electrodes (through nostril, less invasive)
- ☐ Know clinical indications for each electrode type
- ☐ Understand sensitivity differences (sphenoidal > nasopharyngeal > T1/T2 > T3/T4)
- ☐ Know artifact considerations for each electrode type
- ☐ Understand when to use exotic electrodes vs standard system

## 16. Study Tips

1. **Memorize standard 10-20 system first:** Foundation for all electrode placement
2. **Learn exotic electrode locations:** T1/T2, F9/F10, P9/P10 placement techniques
3. **Understand invasiveness:** Sphenoidal (invasive) vs nasopharyngeal (less invasive) vs surface (non-invasive)
4. **Know sensitivity:** Sphenoidal > nasopharyngeal > T1/T2 > T3/T4 for temporal activity
5. **Learn clinical indications:** When to use each electrode type
6. **Understand placement techniques:** Measurement methods for each electrode
7. **Know artifacts:** Common artifacts for each electrode type
8. **ABRET focus:** Expect questions on T1/T2, sphenoidal, nasopharyngeal, and clinical indications

## 17. Internal Cross-Links

### Workflow

- **Electrode Application & 10-20 System:** Standard placement techniques
- **Electrodes & Impedance:** Electrode types and materials
- **Montages & Referencing:** How electrodes are connected

### Patterns

- **Epileptiform Discharges:** Localization using exotic electrodes
- **Temporal Lobe Epilepsy:** When to use T1/T2, sphenoidal, nasopharyngeal

### Cases

- **Temporal lobe epilepsy cases:** Requiring exotic electrodes
- **Presurgical evaluation cases:** Using sphenoidal or nasopharyngeal electrodes

### Quizzes

- **Electrode placement MCQs:** Questions on standard and exotic electrodes
- **Temporal lobe epilepsy questions:** Questions on electrode selection
- **Presurgical evaluation questions:** Questions on invasive electrodes

---

### End of Study Guide

*For additional practice, complete quiz questions tagged: electrode-placement, 10-20, exotic-electrodes, T1-T2, sphenoidal, nasopharyngeal, temporal-epilepsy*